## Chemical Indicators

## Results Reference Guide

Unprocessed

Failures

One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again. The main failure causes are:

- Inadequate cycle duration
- Inadequate cycle selection
- Overloading / incorrect loading of the sterilizer
- Improper sterilization agent concentration
- Sterilizer malfunction
- Blocking indicator with large metal mass

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## Chemical Indicators CD40

For Hydrogen Peroxide sterilization processes

Type 4 Indicators according to ISO 11140-1:2014

## Product description

Chemdye® CD40 chemical indicators have been designed to react to Plasma or Vaporized Hydrogen Peroxide sterilization processes within loads, ensuring an adequate control of the effectiveness of sterilization processes. Purple ink was developed to change to green when the process reaches the stated values of a critical process variable.

## Characteristics

Chemical strips (105 x 18 mm) printed with indicator ink.  
Indicators 100% free of toxic heavy metals.

## Indications for use

Use Chemdye® CD40 indicators for monitoring Plasma or Vaporized Hydrogen Peroxide sterilization processes.  
Conditions: 6 minutes, 50 °C, 2,3 mg/l H2O2.

## Precautions

Do not expose Chemdye® CD40 indicators to Steam, Ethylene Oxide or any sterilization process other than Plasma or Vaporized Hydrogen Peroxide.

## Instructions for use

- Place a Chemdye® CD40 indicator within each package, bag or tray that will be sterilized by Plasma or Vaporized Hydrogen Peroxide.
- In rigid containers place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.
- After completing the process, analyze results. The indicator must turn to the reference color for considering that the indicated conditions were met.

## Storage and shelf life

- Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.
- Indicators have an expiration date of 5 years from the date of manufacture when stored under recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.
- After use, the indicator will not change within 1 year when stored at previously indicated conditions.

## Disposal

Discard chemical indicators after use according to your country's healthcare and safety regulations