Multi-variable Indicators

Results Reference Guide

Unprocessed

Failures

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Multi-variable Indicators CD50

For Formaldehyde sterilization processes

EN

Type 4 Indicators according to ISO 11140-1:2014

Product description

Chemdye® multi-variable indicators are designed to react to Low Temperature Steam and Formaldehyde (LTSF) sterilization processes within loads, ensuring an adequate control of the effectiveness of sterilization processes (temperature, time and gas concentration). Bordeaux/red ink was developed to change to green when the process reaches the stated values of a critical process variable. Multi-variable indicators allow to rapidly detect whether the critical parameters of the sterilization processes for which they have been designed to respond, have been reached.

Characteristics

CD50: Chemical strips of 105 x 18 mm.
Indicators 100% free of toxic heavy metals.

Indications for use

Use Chemdye® multi-variable indicators for monitoring LTSF sterilization processes.

Precautions

Do not expose Chemdye® multi-variable indicators to Steam, Ethylene Oxide or any sterilization processes other than LTSF.

Instructions for use

- Place a Chemdye® multi-variable indicator within each package, bag or tray that will be sterilized by LTSF.
- In rigid containers place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.
- After completing the process, analyze results. The indicator must turn to the reference color for considering that the indicated conditions were met.

Storage and shelf life

- Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.
- Indicators have an expiration date of 5 years from the date of manufacture when stored under recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.
- After use, the indicator will not change within 1 year when stored at previously indicated conditions.

Innovation in technologies for sterilization and disinfection control

Industria Argentina - Made in Argentina
Fabricado por Terragene S.A. - Güemes 2879 - (2000)
Rosario - Santa Fe - Argentina

Multi-variable Indicators

Results Reference Guide

Unprocessed

Failures

One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again.

The main failure causes are:

- Inadequate cycle duration
- Inadequate cycle selection
- Overloading / incorrect loading of the sterilizer

- Improper sterilization agent concentration
- Sterilizer malfunction
- Blocking indicator with large metal mass

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