Unique Point Integrator Indicator


Quality Certification

STEAM

Product code _____

Temperature _____ SV (Time) _____
 121 °C _____
 128 °C _____
 135 °C _____

SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

According to ISO 11140-1:2014 standard

Lic. Adrián J. Rovetto
 Director Técnico
 Quality Assurance Director

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Type 5 Indicator according to ISO 11140-1:2014

Product description

INTEGRON® unique point integrator indicators were developed for monitoring Steam sterilization processes between 121 °C and 135 °C, ensuring an adequate control of the effectiveness of sterilization processes (temperature, time, quality of steam). Yellow ink was developed to turn to black when a theoretical spores population reaches its kill time, indicating integration condition has been reached. This condition is calibrated with the kill time of a 10⁶ *Geobacillus stearothermophilus* spores population, calculated in BIER (Biological Indicator Evaluator Resistometer).

Characteristics

IT26-1YL and IT26-1YS consist of paper strips, 142 mm x 23 mm and 72 mm x 23 mm in size, respectively.

Indications for use

Use the integrator indicator to monitor Steam sterilization processes between 121 °C and 135 °C. Integration condition is reached in Steam sterilization processes as described in the quality certificate.

Precautions

Do not expose this product to EO, Dry Heat, or any sterilization process other than Steam.

Directions for use

- Place an INTEGRON® unique point Integrator Indicator within each package, bag or tray that will be sterilized by Steam.
- In rigid containers place an INTEGRON® unique point integrator indicator in every corner, or at least two diagonally opposite corners, to ensure that steam has reached the most inaccessible areas.
- After completing the sterilization process, analyze the results. If the ink printed on the indicator changes to the reference color printed on the product, integration condition is met and the sterilization process was successful.

Storage and shelf life

- Indicators should be stored in dry places, away from light, at temperatures between 10-30 °C and 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.
- Indicators have an expiration date of 5 years from the date of manufacture when stored under the recommended conditions. The expiration date and the batch number of the product are specified on the package label and in each unit of product. Do not use indicators after their expiration date.
- After use, the indicator will not change for 1 year if stored at temperature and relative humidity previously indicated.

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


Failures

Processed and correct

One or more parameters of sterilization were not achieved. The materials cannot be considered sterile and must be processed again. These results are usually obtained when steam quality is poor or adequate vacuum has not been generated.

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