

Indicador emulador de un punto Unique point emulator indicator

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Results reference guide

	↔		Unprocessed
	↔		Failures One or more parameters of sterilization were not achieved. The materials cannot be considered sterile and must be processed again. These results are usually obtained when steam quality is poor or adequate vacuum has not been generated.
	↔		
	↔		
	↔		Processed and correct

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Class 6 indicator according to EN ISO 11140-1:2009

Product description

The Integron® indicator emulators are designed to react in steam sterilization processes inside of the loads, ensuring an adequate control of the effectiveness of sterilization processes (steam, time, quality of steam). The yellow ink was developed to change to black when the process reaches the variables set out in the emulator.

Characteristics

Emulator indicators IT27-3YS, IT27-4YS, IT27-5YS, IT27-7YS, IT27-9YS and IT27-18YS consist of a paper strip of 72 x 23 mm.

Emulator indicators IT27-3YL, IT27-4YL, IT27-5YL, IT27-7YL and IT27-9YL consist of a paper strip of 142 mm x 23 mm.

Indications for use

Use Integron® indicator emulators for monitoring Steam sterilization processes.

Emulation condition:

IT27-3YS	3 min. 134 °C
IT27-3YL	3 min. 134 °C
IT27-4YS	4 min. 134 °C
IT27-4YL	4 min. 134 °C
IT27-5YS	5 min. 134 °C - 15 min. 121 °C
IT27-5YL	5 min. 134 °C - 15 min. 121 °C
IT27-7YS	7 min. 134 °C - 20 min. 121 °C
IT27-7YL	7 min. 134 °C - 20 min. 121 °C
IT27-9YS	9 min. 134 °C
IT27-9YL	9 min. 134 °C
IT27-18YS	18 min. 134 °C

Precautions

Do not expose these products to EO, dry heat, or any sterilization processes other than Steam. Do not reuse the sterilizer until the emulator's color reaches the emulation condition.

Instructions of use

- Place an Integron® indicator emulator within each package, bag or tray that will be sterilized by steam.
- In rigid containers place an indicator emulator in every corner, or at least two diagonally or opposite corners, to ensure that the steam has reached the most inaccessible areas.
- After completing the process, analyze the results. The indicator must change to the reference color for considering that emulation conditions were reached.

Storage and shelf life

- Indicators should be stored in dry places, protect from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not storage close to sterilizing agents.
- Indicators have an expiration date of 5 years from the date of manufacture when stored in the recommended conditions. Both expiration date and batch of the product are specified on the package label and on each unit of product. Do not use indicators after their expiration date.
- After use, the indicator will not change in 1 year if stored at temperature and relative humidity previously indicated.

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